

LC-Resin™

Easy & Fast, Light cured pattern resin

One-Step, Connection, Perfect Fit



No mixing. No odor. No shrinkage.

The ideal solution for precise and efficient prosthetic connection.

LC Pattern Resin

Light-Cured Precision for Perfect Connection.



MEDICLUS CO., LTD.

LC-Resin™

Save Time, No Smell, No Shrinkage, No Ash

LC-Resin™ offers a convenient syringe-type design for quick and effortless application. Its one-step, pre-mixed formulation eliminates the need for powder-liquid mixing, allowing for faster, cleaner procedures. The material is odor-free, ensuring a pleasant working environment, and leaves no residue upon removal, minimizing material waste. With zero shrinkage (0%), it provides outstanding dimensional stability and accuracy. In addition, LC-Resin™ is highly cost-effective, making it an ideal solution for a wide range of clinical cases.

Advantages

- Fast working time with light curing type
- Flexural strength : 50Mpa
- Elastic modulus : 1200Mpa
- Curing depth : 5mm


Usage

- Build a jig between implant impression copings
- Build a pattern for bridge, inlay and onlay
- Build a pattern for metal framework of denture
- Build a pattern for clasp of partial denture
- Build a outer shell of telescope

Application



Implant jig






Gold crown



Cap crown

Application process



Delivery Forms



LC-Resin™

#602-1 LC-Resin Red 14g x 1EA / Disposable tips

#602-2 LC-Resin Violet 14g x 1EA / Disposable tips

MEDICLUS

Headquarter | No.1210, 134, Gongdan-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28576, Korea | +82 43 211 2877

Sales Office | No.206 Brown Stone, 464, Cheongpa-ro, Jung-gu, Seoul, 04510, Korea | +82 10 8385 2877

www.mediclus.co.kr | sales@mdclus.com

