

Any-Paste®

Water-Soluble Calcium Hydroxide Temporary Root Canal Filling Material

Easy Application Convenient Removal



Versatile Application

Applicable throughout all stages of endodontic treatment and for a wide range of periapical conditions.

Multifunctional clinical benefits — antibacterial activity, exudate control, and more.

Biological Stabilization Support

Prevents bacterial contamination and inhibits microbial proliferation during endodontic treatment.

Dissolves residual pulp tissue within the canal to help relieve pain.

MDCLUS
MEDICLUS CO., LTD.

Any-Paste®

Essential Calcium Hydroxide Temporary Root Canal Filler for Endodontic Treatment

A pre-mixed calcium hydroxide paste for temporary root canal filling that helps control infection, prevent bacterial regrowth, and ensure a clean, stable canal environment. Promotes healing in exudative canals and provides effective pain relief by dissolving remaining pulp tissue.

Any-Paste® features excellent flowability and a clinically friendly formulation, making it useful for infection control and evaluation of tissue condition prior to the application of One-Fil.


Indications

- Apexification procedures
- Temporary root canal filling
- Direct pulp capping and vital pulp therapy
- Pulpotomy
- Management of persistent intracanal exudation
- Management of internal and external root resorption

Key Features & Benefits

- Convenient and precise application — ready-to-use premixed formulation for quick and accurate placement.
- Excellent flowability — ensures smooth and uniform intracanal application.
- Outstanding antibacterial effect — high alkalinity (pH 12.6) effectively inhibits bacterial growth.
- High radiopacity (400% Al) — provides excellent radiographic visibility.
- Easy removal — hydrophilic formulation allows effortless removal with EDTA or NaOCl irrigation.

Clinical Images




Before treatment



Application of
Any-Paste®



After removal



After One-Fil®
application

Delivery Forms



#505-1 Any-Paste®

Any-Paste® 2g x 1ea

MEDICLUS

Headquarter | No.1210, 134, Gongdan-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28576, Korea | +82 43 211 2877

Sales Office | No.206 Brown Stone, 464, Cheongpa-ro, Jung-gu, Seoul, 04510, Korea | +82 10 8385 2877

www.mediclus.co.kr | sales@mdclus.com

